DK50 2x4VR/110/M

- EN User manual
- DE Benutzerhandbuch
- FR Manuel de l'utilisateur
- RU Руководство пользователя
- PL Instrukcja obsługi
- SK Návod na použitie
- CS Návod k použití

COMPRESSOR
KOMPRESSOR
COMPRESSEUR
КОМПРЕССОР
SPREŽARKA
KOMPRESOR

DK50 2x4VR/110/M

EKOM spol. s r. o.

Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk

email: ekom@ekom.sk

DATE OF LAST REVISION
DATUM DER LETZTEN ÜBERARBEITUNG
DATE DE LA DERNIÈRE RÉVISION
ДАТА ПОСЛЕДНЕЙ РЕДАКЦИИ
DATA OSTATNIEJ AKTUALIZACJI
DÁTUM POSLEDNEJ REVÍZIE
DATUM POSLEDNÍ REVIZE

07/2021

NP-DK50 2x4VR 110-AD-
3_07-2021
112000447-000

EN	<u>CONTENTS..... 5</u>
DE	<u>INHALT 59</u>
FR	<u>TABLES DES MATIÈRES 115</u>
RU	<u>ОГЛАВЛЕНИЕ 171</u>
PL	<u>SPIS TREŚCI 229</u>
SK	<u>OBSAH 283</u>
CS	<u>OBSAH 337</u>

CONTENTS

IMPORTANT INFORMATION	6
1. CONFORMITY WITH THE REQUIREMENTS OF THE EUROPEAN UNION.....	6
2. INTENDED USE	6
3. CONTRAINDICATIONS AND SIDE-EFFECTS	6
4. SYMBOLS.....	6
5. WARNINGS	7
6. STORAGE AND TRANSPORT CINDITIONS.....	9
PRODUCT DESCRIPTION	10
7. VARIANTS	10
8. ACCESSORIES	11
9. PRODUCT FUNCTION	11
TECHNICAL DATA.....	16
INSTALLATION.....	21
10. INSTALLATION CONDITIONS	21
11. COMPRESSOR ASSEMBLY	21
12. PNEUMATIC CONNECTION	24
13. ELECTRICAL CONNECTION	28
14. COMMISSIONING	32
15. PNEUMATIC AND ELECTRICAL DIAGRAMS	33
OPERATION.....	36
16. SWITCHING ON THE COMPRESSOR.....	37
17. SWITCHING OFF THE COMPRESSOR.....	37
18. AD DRYER OPERATING MODES:.....	37
PRODUCT MAINTENANCE	40
19. PRODUCT MAINTENANCE.....	40
20. LONG-TERM SHUTDOWN.....	55
21. DISPOSAL OF DEVICE	55
TROUBLESHOOTING.....	56
22. REPAIR SERVICE	58
ANNEX	391
23. INSTALLATION RECORD	391

IMPORTANT INFORMATION

1. CONFORMITY WITH THE REQUIREMENTS OF THE EUROPEAN UNION

This product conforms to the requirements of the Regulation (EU) on medical devices (MDR 2017/745) and is safe for the intended use if all safety instructions are followed.

2. INTENDED USE

The compressor is used as a source of clean, oil-free compressed air to power active medical devices where the parameters and properties of the compressed air are suitable for the specific application.

Compressed air supplied by the compressor is unsuitable for use with artificial lung ventilation devices without further filtration

Any other use of the product beyond its intended use is considered an incorrect use. The manufacturer is not liable for any damages or injuries resulting from the incorrect use.

3. CONTRAINDICATIONS AND SIDE-EFFECTS

There are no contraindications or side-effects known.

4. SYMBOLS

The following symbols and marks are used in the User manual, on the device and its packaging:

General warning

Warning - risk of electric shock

Warning - compressor is controlled automatically

Warning - hot surface

General caution

Read the operating instructions

Refer to instruction manual

CE – marking

Serial number

Protecting earthing

Fuse

	Compressed air inlet
	Compressed air outlet
	Control wire input
	Package handling label – fragile
	Package handling label – this side up
	Package handling label – keep dry
	Package handling label – temperature limits
	Package handling label – limited stacking
	Package label – recyclable material
	Manufacturer

5. WARNINGS

The product is designed and manufactured to be safe for the user and the surrounding environment when used in the defined manner. Keep the following warnings in mind.

5.1. General warnings

READ THE USER MANUAL CAREFULLY BEFORE USING THE DEVICE AND KEEP IT FOR FUTURE USE!

- The user manual supports the correct installation, operation and maintenance of the product. Careful review of this manual will provide the information necessary for the correct operation of the product for its intended use.
- Original packaging should be retained for the return of the device. Only the original packaging ensures protection of the device during transport. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.
- This warranty does not cover damages originating from the use of accessories other than those specified or recommended by the manufacturer.
- The manufacturer only guarantees the safety, reliability and function of the device if:
 - installation, new settings, changes, expansion, and repairs are performed by the manufacturer or an organisation authorised by the manufacturer.
 - the product is used pursuant to the user manual.
- The user manual corresponds to the configuration of the product and its compliance with the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights for the protection of its configuration, methods and names.

- Translation of the user manual is performed in accordance with the best available knowledge. The Slovak version is to be used in the event of any uncertainties.
- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

5.2. General safety warnings

The manufacturer designed and manufactured the product to minimise all risks when used correctly for the intended use. The manufacturer considers it its obligation to lay down the following general safety precautions.

- Use and operation of the product must comply with all laws and local regulations valid in the place of use. The operator and user are responsible for following all the appropriate regulations in the interests of performing work safely.
- Only the use of original parts guarantees the safety of operating personnel and the flawless operation of the product itself. Only the accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer should be used.
- The operator must ensure that the device is functioning correctly and safely before every use.
- The user must be familiar with the operation of the device.
- Do not use the product in environments with a risk of explosion.
- The user must inform the supplier immediately if any problem directly related to the operation of the device occurs.

5.3. Safety warnings on protection from electric current

- The device must only be connected to a properly installed, earthed socket.
- Before the product is plugged in, ensure that the mains voltage and frequency stated on the product are in compliance with the values of the mains.
- Prior to putting the device into operation check for any damage to the connected pneumatic lines and electrical wiring. Replace damaged pneumatic lines and electrical wirings immediately.
- Immediately disconnect the product from the mains (remove power cord from the socket) in hazardous situations or in the case of a technical malfunction.
- During all repairs and maintenance, ensure that:
 - the mains plug is removed from the power socket
 - pressure is vented from the air tank and pipes
- The product shall only be installed by a qualified technician.

6. STORAGE AND TRANSPORT CINDITIONS

The compressor is shipped from the manufacturer in transport packaging. This protects the product from damage during transport.

Potential for damage to pneumatic components.

The compressor must be transported only when all air has been vented. Before moving or transporting the compressor, release all the air pressure from the tank and pressure hoses, from dryer chambers and drain condensate from the tank and from the condensate separator on the dryer.

Keep the original factory packaging in case the device needs to be returned Use the original factory packaging during transport as it provides optimum protection for the product. . If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.

The compressor is shipped in a vertical position and must be secured using transport straps.

Protect the compressor from humid and dirty environments and extreme temperatures during transport and storage. Do not store near any volatile chemical substances.

If not, please dispose of the original packaging material in an environmentally-friendly way. The packaging cardboard can be recycled with old paper.

Storing or shipping the equipment in any conditions other than those specified below is prohibited.

6.1. Ambient conditions

Products may only be stored and transported in vehicles that are free of any traces of volatile chemicals under the following climactic conditions:

Temperature

–25°C to +55°C, 24 h at up to +70°C

Relative humidity


max. 90% (non-condensing)

PRODUCT DESCRIPTION


7. VARIANTS

The compressor is manufactured according to its intended application in the following variants:

- DK50 2x4VR/110/M** Compressor with adsorption dryer
DK50 2x4VR/110S/M Compressor in cabinet with adsorption dryer
- Cabinet S110R** It decreases noise level of the compressor.

DK50 2x4VR/110/M

S110R
DK50 2x4VR/110S/M

8. ACCESSORIES

Accessories that are not included in the standard order must be ordered separately.

Set of compressed air outlet filters

The compressor may be equipped with a set of filters if specified. The filter set may be equipped with an air pressure regulator.

Where a different level of air filtration is required, this requirement must be agreed upon with the supplier and specified in the order.

Type	Use	Level of filtration (µm)	Bypass function *	Article number
FS 40F		1		604014119-000
FS 40M	DK50	1+0.1	no	604014119-004
FS 40S	2x4VR/110	1+0.01		604014119-024
FS 40AH		1+AC+HC (0.01)		604014119-005

*) These FS do not contain a filter bypass, which will ensure a continuous flow of air when replacing the filter element. Such a set must be ordered separately.

Filter set regulator assembly

The compressor may be equipped with a pressure regulator of the compressed air outlet if specified. The regulator must be selected according to the application to the filter set, or separately. The regulator shall ensure constant pressure at the outlet.

Type	Use	Article number
Regulator complete		604014125-000
REG15	DK50 2x4VR/110/M	447000001-243

Filter set brackets

A suitable bracket must be ordered for every filter set.

Type	Use	Article number
Compressor-mounted bracket		604014131-000
Wall-mounted bracket	DK50 2x4VR/110/M	603014120-000

9. PRODUCT FUNCTION

9.1. Compressor with adsorption dryer

The compressor aggregate (1) draws in air through the inlet filter (8) and compresses it through a check valve (3) and into a manifold, from which it is routed to the adsorption dryer (9) through a connecting hose. From the inlet to the dryer module, the air is first cooled in the integrated cooler (14) and then moves through the condensate separator (40), entering the active chamber with adsorbent (16), where the air is dried. A portion of the air is fed into the second, regeneration chamber, where this air is used to remove moisture from the adsorbent, after which it is released through the silencer (38). The active chamber switches on a cyclical basis. The dry and filtered air then passes through the check valve (3) into the air tank (2).

9.2. Compressor cabinet

The soundproof cabinet provides compact covering of the compressor, dampens noise efficiently, providing sufficient exchange of cooling air. The fan (10) under the compressor pump provides cooling of the compressor. It is in operation at the same time with the compressor motor or after switching on the temperature switch at temperature above 40°C. After cooling the cabinet area approximately under 32°C, the fan turns off automatically.

Risk of compressor overheating.

Make sure that there are no obstacles at the cooling air inlet into the cabinet (around the bottom part of the cabinet) and at the hot air outlet on the top back side of the cabinet.

If the compressor is placed on a soft floor, e.g. carpet, create space between the base and the floor or the cabinet and the floor, e.g. underlay the footings with hard pads to ensure sufficient cooling of the compressor.

Descriptions for figures 1 - 3

1	Air pump	23	Elbow fitting
2	Air tank	24	–
3	Non-return valve	25	Compressor wheel
4	Pressure switch	26	Cabinet fan
5	Safety valve	27	Handle
6	Pressure gauge	28	Switch
7	Drain valve	29	Maintenance indicator
8	Inlet filter	30	Cabinet pressure gauge
9	Adsorption dryer	31	-
10	Compressor fan	32	-
11	Magnetic holder	33	Regeneration solenoid valve
12	Condensate collection vessel	34	Regeneration solenoid valve
13	Compressed air outlet	35	Inlet solenoid valve
14	Cooling module	36	Inlet solenoid valve
15	Inlet valve module	37	Relief valve
16	Dryer chamber	38	Noise silencer
17	Outlet module	39	Compressed air inlet
18	Control unit	40	Condensate separator
19	Cabinet gas springs	41	Pressure gauge
20	Wall stop	42	Automatic condensate drain
21	Power supply cord	43	Dryer pan
22	Pressure gauge hose	44	Air outlet

Fig. 1: DK50 2x4VR/110/M – Compressor with adsorption dryer


Fig. 2: Adsorption dryer




Fig. 3: Cabinet

TECHNICAL DATA

Compressors are designed for operation in dry, ventilated and dust-free indoor rooms under the following climactic conditions:

Temperature +5°C to +40°C

Relative humidity max. 70%

Working pressure 6 – 8 bar		DK50 2x4VR/110/M	DK50 2x4VR/110S/M
Nominal voltage	V, Hz	3x400V, 50Hz	3x400V, 50Hz
Frequency ^{a)}			
Capacity at 6 bar (FAD) -20°C	l/min	450	450
Working pressure ^{b)}	bar	6.0 – 8.0	6.0 – 8.0
Max. current	A	9.0	9.8
Motor power	kW	2x2.2	2x2.2
Air tank volume	l	110	110
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	9.0	9.0
Noise level at 5 bar (L _{pA})	dB	≤76.7	≤58.0
Operating mode	%	S1-100	S1-100
PDP drying performance at 7 bar ^{d)}	°C	≤-20	≤-20
Time to fill air tank from 0 to 7 bar	s	88	88
Dimensions (net) w x d x h	mm	2000x503x965	2000x834x1045
Net weight ^{c)}	kg	186 ^{c)}	268 ^{c)}
Classification under EN 60601-1	Class I.		


^{a)} Specify the compressor version when ordering

^{b)} Consult any other range of pressure with the supplier

^{c)} Weight is indicative and only applies to the product without accessories

^{d)} Applies to ambient temperatures of <30°C PDP – pressure dew point


Dependence of compressor output on working pressure

Working pressure 8 – 10 bar		DK50 2x4VR/110/M	DK50 2x4VR/110S/M
Nominal voltage Frequency ^{a)}	V, Hz	3x400, 50	3x400 , 50
Capacity at 8 bar (FAD) -20°C	l/min	345	345
Working pressure ^{b)}	bar	8.0 – 10.0	8.0 – 10.0
Max. current	A	9.4	10.2
Motor power	kW	2x2.2	2x2.2
Air tank capacity	l	110	110
Air quality - filtration	µm	-	-
Maximum operating pressure of safety valve	bar	12.0	12.0
Noise at 5 bar (L _{pA})	dB	≤76.7	≤58.0
Operating mode	%	S1-100	S1-100
PDP drying performance at 7 bar ^{d)}	°C	≤ -20	≤ -20
Time to fill air tank from 0 to 7 bar	s	121	121
Dimensions (net) w x d x h	mm	2000x503x965	2000x834x1045
Net weight ^{c)}	kg	186 ^{c)}	268 ^{c)}
Classification under EN 60601-1		Class I.	

- ^{a)} Specify the compressor version when ordering
- ^{b)} Consult any other range of pressure with the supplier
- ^{c)} Weight is indicative and only applies to the product without accessories
- ^{d)} Applies to ambient temperatures of <30°C PDP – pressure dew point

Dependence of compressor output on working pressure

FAD correction of capacity for altitude

Capacity given in the form of FAD („Free Air Delivery“) applies to the following conditions:

Altitude	0 m.n.m.	Temperature	20°C
Atmospheric pressure	101325 Pa	Relative humidity	0%

To calculate FAD compressor capacity in dependence on altitude, it is necessary to apply correction factor according to the following table:

Altitude [m.n.m.]	0 -1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD correction factor	1	0.80	0.71	0.60

Electromagnetic compatibility declaration

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

Guidance and manufacturer's declaration - electromagnetic emissions		
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	$U_T=0\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=0\%$, 1 cycle $U_T=70\%$ 25/30 cycles (at 0°) $U_T=0\%$, 250/300 cycles	$U_T=>95\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=>95\%$, 1 cycle $U_T=70\%$ (30% dipp U_T), 25(50Hz)/30(60Hz) cycles (at 0°) $U_T=>95\%$, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.
Power frequency (50/60 Hz) IEC 61000-4-8	N/A	N/A	The equipment does not contain magnetically sensitive components or circuits.
NOTE: U_T is the A.C. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}, 80 \text{ MHz to } 800 \text{ MHz}$ $d=2,3\sqrt{P}, 800 \text{ MHz to } 2,7 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.			
^b Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.			

INSTALLATION

Risk of incorrect installation.

Only a qualified technician may install the compressor and place it into operation for the first time. Their duty is to train operating personnel on the use and maintenance of the equipment. An entry is made in the equipment installation record to certify installation and operator training. (see the Annex chapter)

10. INSTALLATION CONDITIONS

- The compressor may only be installed and operating in dry, well-ventilated and clean environments under the conditions specified in the Technical Data chapter.

Risk of damage to the device.

The equipment may not be operated outdoors or in otherwise wet or damp environments.

Risk of explosion.

Do not use the equipment in the presence of explosive gases, dust or combustible liquids.

Burn or fire hazard! Caution! Hot surface!

Portions of the compressor, dryer and connecting hoses between the dryer and compressor may be hot and reach hazardous temperatures during compressor operation that may harm materials or operating staff.

- The compressor must be installed so that it is accessible at all times for operating and maintenance. Please ensure that the nameplate on the device is readily accessible.
- The compressor must stand on a flat, sufficiently stable base (be aware of the weight of the compressor, see the Technical Data chapter).

High temperature hazard.

The placement of air flow impediments upstream or downstream of the cooler is prohibited. The temperature of internal and external parts of the cooler may be hot and reach hazardous temperatures.

You may notice a “new product” odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.


11. COMPRESSOR ASSEMBLY

11.1. Handling and releasing the compressor

- Unpack the compressor from the packaging.
- Position the compressor at the site of installation. (Fig. 4)

Use the handles on the compressor to move the compressor as needed. Do not use other parts of the compressor (the aggregate, cooler, etc.) for grip


Fig. 4: Handling the compressor

- Remove the transport stabilisers from the air pumps (Fig. 5).

Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.

Remove all devices used to secure the aggregates once the compressor is installed and levelled at the site of final installation.

Fig. 5: Releasing the air pumps


Assembly of the AD dryer

- Remove the dryer from the packaging.
- Install the dryer in its operating position (Fig. 6).

At least two persons are needed to handle the equipment.

Integrated handles are installed on the lower brackets on the product. Each person must grasp the equipment with one hand on a handle and the other behind the dryer chamber when moving the equipment.

Fig. 6: Handling the dryer

11.2. Placement of the compressor in the cabinet

Opening the upper cover

- Open the lock by turning with a screwdriver as shown in the pictogram and pick up using the handle.
- The gas springs will keep the cover open.

Be careful not to pinch your fingers when closing the cover on the cabinet.

Always lock with the quick releases after closing the cover on the cabinet!

Removal of the front panel on the cabinet

- Remove the door held on by the 4 (four) M5 screws and disconnect the earthing conductor.

Compressor placement

- Insert the compressor into the cabinet so the compressed air outlet is pointing towards the operator.

Cabinet door installation

- Align the door with the cabinet, connect the grounding lead to the door and install the door to the cabinet using 4 screws.

12. PNEUMATIC CONNECTION

12.1. Connecting the dryer

AD dryer compressed air inlet

- Connect the compressed air outlet from the compressor to the dryer inlet (1). Use the shorter of the provided hoses. (2200 mm)
- A G 3/4" connection is installed.


Fig. 7: Compressed air inlet

Fig. 8: Proper hose routing

Burn or fire hazard! Caution! Hot surface!

When installing connecting hose (Fig. 7, Fig. 8) at the air inlet to the dryer, please note the hose and the manifold on the compressor may be hot and reach hazardous temperatures that may harm materials or operating staff.

AD dryer compressed air outlet

- Connect the outlet from the dryer (1) to the air inlet on the air tank. Use the longer of the provided hoses. (2900 mm)
- A G 3/4" connection is installed.


Fig. 9: Compressed air outlet

- When installing the hose, secure the nipple in the dryer with a wrench.

Fig. 10: Compressed air outlet

- Connecting hose (A) is routed from the compressor outlet (1) to the dryer inlet (2) and is routed together with hose B using double clips (3).
- Connecting hose (B) is routed from the dryer outlet (4) to the air tank inlet (5) and is connected to the dryer (6) using clips.


Fig. 11: Routing the connecting hoses

Fig. 12: Location of the first double clip

Condensate outlet from dryer

- Connect a hose to the outlet (1) from the automatic condensate drain (2) to drain piping or to the provided collection vessel.

A noise silencer is recommended when connecting directly to drain piping.


Fig. 13: Condensate drain

12.2. Connecting to the compressed air outlet

- Route the pressure hose from the compressed air outlet (1) to the connected equipment.
- Connect the drain hose (1) to the condensate vessel on compressors with dryers.
- A G 1/2 connection is installed.

Fig. 14: Connecting to the compressed air outlet

12.3. Connecting the cabinet pressure gauge to the compressor

- Remove the plug (1) from the threads (2) on the pneumatic block of the compressor
- Connect the cabinet pressure gauge hose to the threaded fitting


Fig. 15: Connecting the cabinet pressure gauge to the compressor

12.4. Routing hoses and electrical cables in cabinet-mounted compressors

- Route the hose and electrical cables through the opening in the rear wall of the cabinet for compressor models with a dryer.

Risk of damage to pneumatic components.

Air hoses must not be broken.

13. ELECTRICAL CONNECTION

- The product is delivered with a cord equipped with a plug with earthing pin.
- With cabinet-mounted compressors, route the power cord through the opening in the back wall of the cabinet.
- Keep the socket easily accessible to ensure that the device can be safely disconnected from the mains.
- Connection to the power distribution cabinet must be max.16 A

Risk of electric shock.

It is necessary to follow all local electro technical regulations. The mains voltage and frequency must comply with the data stated on the device label.

Risk of fire and electric shock.

Electrical cord must not be broken.

Connecting the control harness and power cord from the compressor

- Connect the control harness and power cord from the compressor to the dryer socket (1).


Fig. 16: Connecting the control harness and power cord

- Connect the pin for equipotential bonding 6Ø mm (1) to the electrical circuit using the defined method under valid electrical engineering regulations
- The socket for equipotential bonding (2) is an accessory and is not provided in the basic package

Fig. 17 :Equipotential bonding

- Mount the connector of the grounding conductor on the door and screw on the cabinet.

13.1. Connecting a compressor not installed in a cabinet

- Insert the mains plug into a rated mains socket.
- The compressor is ready for operation.

13.2. Connecting a compressor installed in a cabinet


- Route the mains plug through the opening in the rear wall of the cabinet for cabinet-mounted compressors.
- Connect the cabinet electrically to the compressor by inserting the provided power cord with connector into an outlet. (Fig. 18)
- Disconnect the cabinet electrically from the compressor by pulling the connector from the outlet while the latch is released.

Fig. 18: Connect the cabinet to the compressor

Ethernet connection

The dryer may be connected to an Ethernet 10/100 M network via the controller as follows:

- Connect the Ethernet network cable to the RJ-45 connector on the rear of the dryer.
- Configuration of the IP address to connect to the local network:
 - The default IP address of the BM module is: 192.168.0.3.
 - The user should request that the service technician change the IP address during installation.

**Fig. 19: Ethernet connection**

Web server

The controller has an integrated Web Server function that facilitates compressor monitoring via a PC, smartphone or tablet using a conventional web browser (Mozilla, Opera, Safari, Google Chrome, etc.)

The process for logging into the Web Server function once the compressor is connected to an Ethernet network is as follows:

- Open the web browser on a PC, smartphone or tablet and enter the IP address of the controller (in this case 192.168.0.3).

- Enter the password “LOGO” and click on the “LOG on” button.

- After logging in, the browser displays the first screen showing the system information for the controller itself: module generation, type, firmware (FW), IP address and activity status.

- Click on the “LOGO! BM” function in the browser to display the current virtual status of the display screen. Navigate through the screen using the ESC and cursor keys the same as on the real display.

Monitoring memory variables

The “LOGO! Variable” function is the other option for monitoring compressor parameters using selected memory variables. Click the “LOGO! Variable” button on the display to bring up a screen and then use the “Add Variable” button to select specific memory variables for monitoring based on the mapping provided by the compressor manufacturer.

Variables assigned based on addresses and models may be viewed in the MAPPING PARAMETERS table (see the Annex chapter).

Select the variable parameters (per the annex) sequentially in the Range, Address, Type and Display Format columns. Variable values are shown in the Value column. The monitoring table may then appear as follows

Variable						
Del Range	Address	Type	Display Format	Value	ModValue	Modify
X VM	0	DWORD	SIGNED	0		✓
X VM	4	WORD	SIGNED	-250		✓
X VM	6	WORD	SIGNED	-500		✓
X VM	8	WORD	SIGNED	-500		✓
X VM	10	DWORD	SIGNED	120000		✓
X VM	14	DWORD	SIGNED	0		✓
X VM	26	WORD	SIGNED	-50		✓
X VM	28	WORD	SIGNED	0		✓

Buttons: Add Variable, Modify All Values

14. COMMISSIONING

- Make sure all transport stabilizers were removed.
- Check that all compressed air hose connections are correct
- Check correct connection to the mains (see chap. 13).
- Check the circuit breaker position, it must be in the "I" position. If the main switch (2) is in the "0" position, turn it to the "I" position. (Fig. 20).
- For cabinet-mounted compressors, turn the main switch (4) on the front side of the cabinet into the "I" position; a green indicator indicates the equipment is in operation. (Fig. 20).
- Check connection of the cabinet connector to the compressor (Fig. 18)
- Check connection of the cabinet manometer hose to the pneumatic block of the compressor (Fig. 15).
- Check to ensure the compressed air connection between the compressor and dryer is correct. (Fig. 11)
- Check to ensure the hose from the automatic condensate drain is connected to the collection vessel.
- Check the connection of the dryer control and power signal wiring from the compressor. (Fig. 16).


The compressor is not equipped with a backup power supply.

15. PNEUMATIC AND ELECTRICAL DIAGRAMS

15.1. Pneumatic diagram

DK50 2x4VR/110/M


Description to pneumatic diagram:


- | | |
|-------------------------|---------------------------------|
| 1 Inlet filter | 14 - |
| 2 Compressor | 15 Relief valve |
| 3 Fan | 16 - |
| 4 OR logic valve | 17 Outlet valve |
| 5 - | 18 Condensate collection vessel |
| 6 Non-return valve | 19 Cooler fan |
| 7 Pressure gauge | 20 Chamber inlet filter |
| 8 Pressure switch | 21 Adsorbent |
| 9 Safety valve | 22 Chamber outlet filter |
| 10 Air tank | 23 Regeneration jet |
| 11 Drain valve | 24 Inlet filter |
| 12 Cooler | 25 Chamber left |
| 13 Dryer solenoid valve | 26 Chamber right |

15.2. Electrical diagrams

DK50 2x4VR/110/M **6 - 8 bar, 8 - 10 bar**
 3/N/PE~400V, 50Hz
 ELEKTRIC MAIN TN-S [TN-C-S]
 ELEKTRIC OBJECT OF 1st. CAT.

DK50 2x4VR/110S/M

Description to electrical diagrams:

M1, M2	Compressor motor	X2, X3	Socket
E1, E2	Compressor fan	B1	Pressure switch
E3, E4	Compressor fan	C10,11	Capacitor
E10, E11	Cabinet fan	F3	Fuse
B2	Temperature switch	P1	Hour counter
Q11, Q12	Contactors	K1, K2	Time relay
		F1,2	Breaker

OPERATION

ONLY TRAINED PERSONNEL MAY OPERATE THE EQUIPMENT!

Risk of electric shock.

In case of emergency, disconnect the compressor from the mains (pull out the mains plug).

Burn or fire hazard.

Portions of the aggregate and compressed air components between the aggregate and the air cooler may be hot and reach hazardous temperatures during compressor operation that may harm materials or operating staff.

Warning – AD dryer is controlled automatically.

Automatic start. The dryer is controlled by a control signal.

Potential for damage to pneumatic components.

The working pressure settings for the pressure switch set by the manufacturer cannot be changed. Compressor operation at a working pressure below the switching pressure indicates high air usage (see the Troubleshooting chapter).

Risk of damage to the dryer.

The dryer may be damaged if operated at ambient temperatures that are higher than the maximum working temperature (ambient temperature $>30^{\circ}\text{C}$ or inlet air temperature of $>100^{\circ}\text{C}$)

Required drying performance can only be achieved when following the defined operating conditions.

Drying performance will decline and the achieved dew point will drop if the dryer is operated at any pressure below the minimum working pressure.

During prolonged operation of the compressor, the ambient temperature around the compressor may increase to over 40°C . At this point the cooling fan automatically switches on. The fan switches off once the space is cooled to around 32°C .

16. SWITCHING ON THE COMPRESSOR

Start the compressor (without a cabinet) at the pressure switch (1) by turning the switch (2) to position "I." This starts the compressor and fills the tank to the switching off pressure, which then shuts off the compressor.

Start the compressor (with a cabinet) at the switch (4) on the front side of the cabinet, the indicator will turn green. This starts the compressor and fills the tank to the switching off pressure, which then shuts off the compressor.

The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air consumption.

Check the values of switch-on and switch-off pressure using the pressure gauge (3). Tolerance of $\pm 10\%$ is acceptable. Pressure in the air tank must not exceed the allowed operating pressure.

Fig. 20: Switching the compressor

Risk of damage to pneumatic parts.

The pressure switch (1) was set up by the manufacturer and any other changes to its settings may be performed only by a qualified technician trained by the manufacturer.

17. SWITCHING OFF THE COMPRESSOR


- Switching off the compressor due to service or any other reason shall be performed using the pressure switch (1) by turning the switch (2) to position „0“ (Fig. 20) and pulling out the mains plug from the socket. This disconnects the compressor from the mains supply.
- Vent the pressure in the air tank to zero by opening the drain valve.
- Open the venting plug (Fig. 25) to vent the pressure in the dryer chambers.

18. AD DRYER OPERATING MODES:

“RUN” mode

The dryer is in “RUN” mode when the control signal from the compressor is active. The cooler fans are activated and the chambers are cyclically switched by the solenoid valves based on the following time diagram T. This is the chamber switching period, and dT represents the pressure balancing time in the chambers prior to switching.

Valve switching schedule – “RUN” mode

"STANDBY" mode

The dryer is in "STANDBY" mode when the control signal from the compressor is inactive. The cooling fans are off and chamber switching is deactivated.


AD500 SE dryer display

Home screen

- The home screen appears for 5 seconds when the main switch S1 on the dryer is switched to the "I" position

"Adsorptiondryer ADxx00Sis ON"

- The display is backlit in white.

The home screen is followed by the RUN MODE and STAND BY MODE screens based on the compressor control signal.

“RUN MODE” screen

- TOTAL HRS – total time the dryer has been energised
- RUN HRS – total dryer cycling time
- TME-to-MT – time until the service interval expires

R U N M O D E														
T	O	T	A	L		H	R	S	0					
						R	U	N		H	R	S	0	
T	I	M	E	-	t	o	-	M	T		5	0	0	0
I	N	L	E	T	-	A								
						P	U	R	G	E	-	B		

“STAND BY MODE” screen

- TOTAL HRS – total time the dryer has been energised
- RUN HRS – total dryer cycling time
- This screen is not backlit

S T A N D B Y												
M O D E												
D R Y E R N O T												
C Y C L I N G												
T	O	T	A	L		H	R	S	0			
R	U	N		H	R	S	0					

Indication of an upcoming service interval

- 100 hours prior to the next service interval, the back lighting changes from white to orange and the display shows the message “SERVICE DUE IN XY HOURS”, where XY indicates the remaining number of hours until service is due.
- RUN HRS – total dryer cycling time (RUN mode)

S E R V I C E D U E														
I	N		5	0	0	0		H	O	U	R	S		
R	U	N		H	R	S	0							

PRODUCT MAINTENANCE

19. PRODUCT MAINTENANCE

The operator should carry out device checks regularly in the intervals defined by applicable regulations. Test results must be recorded.

The equipment has been designed and manufactured to keep maintenance to a minimum. The following work must be performed to preserve the proper and reliable operation of the compressor.

Unauthorised interference hazard.

Repair work outside the framework of standard maintenance (see Chapter 19.1) may only be performed by a qualified technician (an organisation authorized by the manufacturer) or the manufacturer's customer service.

Standard maintenance work (see Chapter 19.1) may only be performed by the operator's trained personnel.

Only use manufacturer-approved replacement parts and accessories.

Danger of injury or equipment damage.

Prior to commencing compressor maintenance, it is necessary to:

- **check if it is possible to disconnect the compressor from the appliance in order to avoid any risk of injury to the person using the appliance or other material damage;**
- **turn off the compressor;**
- **disconnect it from the mains (pulling the cord out of the mains socket);**
- **vent the compressed air from the air tank.**

Prior to commencing dryer maintenance work, first:

- **turn off the compressor and disconnect it from the mains**
- **check the pressure indicator on the dryer and if there is still pressure, it must first be vented from the dryer chamber**

Venting compressed air poses an injury hazard.

Wear eye protection, i.e. goggles, when venting compressed air from the compressed air circuit (air tank) and from the dryer chamber.

Burn hazard.

When the compressor is running or shortly thereafter, certain portions of the air pump, the compressor's compressed air system, parts of the dryer and connecting hoses between the compressor and dryer may be hot - do not touch these components.

The removed grounding conductor during service must be connected back to the original position after completing the service.

- For any maintenance or repair work, the compressor may be removed from the cabinet by pulling it from the cabinet (on rolling casters as much as the cord between the compressor and the cabinet allows). (see chap. 19) The required maintenance or repair work may then be performed.
- If the cord between the compressor and the cabinet becomes disconnected (by pulling the mains plug, see chap. 13.2), the compressor shall not be functional. In order to bring the compressor into function again, please follow, see chap. 13.2 or 19.8 Or 19.9.

19.1. Maintenance intervals

Time interval	50 Hz	Once a day	Once a week	Once a year	Once every 2 years	2000 h	4000 h	6000 h	8000 h	10000 h	12000 h	20000 h	Chap.	Set of replacement parts	Perform ed by	
															operator	qualified technician
Check of product operation		x											19.2	-		
Drain condensate from air tank – at high humidity		x											19.5	-		
Drain condensate from the air tank - at normal humidity			x										19.5	-		
Check of product function			x										9	-		
Check of pneumatic connections leakage and device inspection				x									19.3	-		
Inspection of electrical connections				x									19.4	-		
Check of cooler and fan				x									19.16	-		
Replacement of the dryer's internal filters										x			19.12	025200322-000		
Check of safety valve				x									19.6	-		
Replacement of pump inlet filter and pre-filter 2V ^{a)}					x			x	x		x		19.7	025200139-000 025200150-000		

^{a)} data is in hours, if not available, data is in years

^{b)} check the operation of the relief valve on the dryer after every dryer malfunction

Time interval	50 Hz						Once every 2 years		Chap.	Set of replacement parts	Perform ed by
	Once a day	Once a week	Once a year	2000 h	4000 h	6000 h	8000 h	10000 h			
Replacement of cassettes with adsorbent media AD500 SE dryer									19.13	603031810-000	qualified technician
Replacement of the dryer's logic valve ball									19.14	074000254-000	
Replacement of the dryer's silencer								x	19.15	025400339-000	
Replacement of the dryer's NC solenoid valve									19.17	025300117-001	

a) data is in hours, if not available, data is in years

b) check the operation of the relief valve on the dryer after every dryer malfunction

19.2. Check of product operation

- Check aggregate condition – the aggregates should be operating normally without excessive vibration or noise. Troubleshoot any problem or call in service personnel if trouble is detected.
- Visually inspect fan operation – the fans must be operating when the aggregates are running. Troubleshoot any problem or call in service personnel if trouble is detected.
- Check to ensure the power cord, the connecting compressed air hoses are undamaged. Replace damaged components or call in service personnel.
- Check the ambient temperature – the ambient temperature must be below the temperature limit (40°C). Cool the space if the temperature is high.
- Compressor with air dryer - open the plug on the condensate collecting vessel and drain the condensate.
- Check the service indicator on the display screen and schedule service as needed.
- Check the operating condition of the equipment (see Chapter 19.4).

19.3. Check the compressed air connections for leaks and inspect the equipment

Leak testing

- Check the compressor's compressed air lines for leaks during operation – pressure supplied by the compressor.
- Use a leak analyser or soapy water to check all joints and connections for leaks. Tighten or reseal the connection where leaks are found.

Inspecting the equipment

- Check the condition of the compressor aggregate for normal operation and noise levels.
- Fan operation check – the fans must be running during the defined compressor work cycles.
- Check function of the temperature switch (B2) – heat the temperature switch to temperature above 40°C (e.g. using a heat gun – Be careful, do not operate with high temperature on plastic components in the surrounding, it can cause deformation of plastic materials). Fan E1, E3 - the compressor must be under voltage.
- Check filters condition – filters must be undamaged and sufficiently clean.
- Check condition of the pump itself, check if there is no dirt inside the crankcase or clearance in the crankshaft.
- Replace any defective parts if there are any failures.
- Replace any defective parts as needed.
- Check the solenoid valves in the valve module – the valves should cyclically cycle between the chambers based on the description of the operation.
- Check the operating condition of the equipment (see Chapter 19.4).

19.4. Inspection of electrical connections

Risk of electric shock.

Inspect the product's electrical connections when the mains are disconnected.

- Check mechanical function of the main switch.
- Check if the power cable and conductors are not damaged.

- Visually check if cables are connected to the terminal box.
- Check all screw connections of the green-yellow PE grounding conductor.

Checking the operating condition of the AD dryer

- Check the service indicator and alarms on the display screen as specified in Chapter 18 and schedule service or repairs as needed.
- Only a service technician is permitted to modify program parameters using the cursor buttons on the LOGO! module and only after a password is entered.
- The operator may monitor the individual values on the display but has no access to the parameters of the program.

Indication of an upcoming or expired service interval

Configuration of a new service interval

Press ESC + ► together and hold for 10 seconds to reset the service interval. Once complete, the display returns back to the home screen.

The SET UP function for a new service interval is only functional if the dryer indicates an upcoming service interval or that a service interval has expired. Once complete, the MAINT counter increases (number of completed service operations) by one. See the screen: Counters.

Information screens – version AD500 S

Information screens are activated by simultaneously pressing and holding down the ESC key and cursor buttons.

- ESC+◀- “COLUMNS TIMING” screen

- Information about the defined valve cycling times
- INLET_A(INLET_B) – duration of the adsorption phase (half-cycle) for chamber A (chamber B)
- PURGE_A(PURGE_B) – duration of the regeneration phase for chamber A (chamber B)

C O L U M N S	
T I M I N G	
P U R G E _ A =	1 0 0 s
I N L E T _ A =	1 2 0 s
I N L E T _ B =	1 2 0 s
P U R G E _ B =	1 0 0 s

- ESC+▼ – “COUNTERS” screen

- Information about counters, the number of times service has been performed and the number of times the solenoid valves have been activated
- MAINT – number of service operations
- VALVE – number of times the solenoid valves have been activated (the number displayed indicates the number of times the Inlet_A inlet valve into chamber A has been activated)
- Version – software version

C O U N T E R S :	
M A I N T	0 x
V A L V E	4 x
V e r s i o n	V 3 . 0 0

- ESC+▲ – “OPERATING TIME” screen
 - information on operating hours: TOTAL, RUN, STAND BY, time remaining to the next service interval and value of the defined service interval
 - TOTAL HRS – total time the dryer has been energised
 - RUN HRS – total dryer cycling time
 - STAND.HRS – total time in STAND BY mode
 - TIME-to-MT – time remaining to the next service interval
 - SERV.INTER – value of the defined service interval

T	O	T	A	L	H	R	S		0
									0
S	T	A	N	D	.	H	R	S	0
T	I	M	E	-	t	o	-	M	T
									1 h
S	E	R	V	.	I	N	T	E	R
									0 h
S	E	T	U	P	=	E	S	C	+ ▶ / 1 0 s /

- ESC+▶ (press and hold for 10 seconds) – to configure a new service interval once a service operation is complete (see above – Configuration of a new service interval).

19.5. Condensate drain

Open the drain valve carefully and slowly. If the drain valve is opened too much, it is accompanied by high noise level and uncontrolled jet of the accumulated condensate.

A wet floor resulting from overflow from the vessel poses a slip hazard.

Condensate from compressors with air dryers is automatically drained into a vessel to collect condensate.

- Monitor the level in the vessel using the markings (depending on the volume of the vessel), and empty at least once a day.


Fig. 21: Check of the condensate collection vessel

Before the following checks it is required:

- For compressor variant with cabinet – unlock the door lock and open the cabinet door.

19.6. Check of safety valve

- Turn the screw (2) on the safety valve (1) several times to the left until the safety valve releases the air.
- Let the safety valve vent for a few seconds.
- Turn the screw (2) fully to the right, the valve must be closed now.

Fig. 22: Check of safety valve

Damage to the safety valve could cause pressure to rise to hazardous levels.

Never use the safety valve to release the air pressure in the air tank. This could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer.

Never adjust a safety valve.

Venting compressed air poses an injury hazard.

Wear safety glasses when inspecting a safety valve.

19.7. Inlet filter replacement

Inlet filter replacement:

- Remove the locking spring (3).
- Pull out the rubber plug by hand (2).
- Remove the dirty intake filter (1).
- Insert a new filter and replace the rubber plug.
- Secure the plug with the locking spring (3).

Pre-filter replacement:

- Pull out the pre-filter by hand (4).
- Replace with a new one and insert it back.


Fig. 23: Inlet filter replacement

19.8. Procedure for connecting a compressor disconnected from the cabinet

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

The compressor (without cabinet) requires that the jumper is always mounted in the terminal strip for correct operation (Fig. 24 pos. D). It then replaces the switch on the cabinet.

If there is no jumper in the terminal strip of the compressor (disconnected from the cabinet), the compressor will not work!


If the cord between the compressor and the cabinet becomes disconnected (by pulling the mains plug) and the compressor is removed from the cabinet, the compressor will not be functional. This is why it is necessary to first create a connection to the terminal strip with a jumper (replacing the circuit breaker function) (Fig. 24 pos. A).


PROCESS:

Jumper installation (A-B-C-D process)

- Disconnect the product from the mains by pulling the plug out of the socket.
- Remove the cover on the electrical panel (on the compressor).
- The jumper is not in the terminal strip – A.
- Lift up the fuse cover – B.
- Insert the jumper into the terminal strip – C.
- Close the fuse cover – D.
- Reinstall the cover on the electrical panel.
- Connect the device to the electrical mains.
- Activate the compressor by turning on the switch on the pressure switch.

Fig. 24: DK50 2x4VR/110

400 V

19.9. Procedure for connecting a compressor to a new cabinet

Prior to any maintenance or repair work, switch off the compressor and disconnect it from the mains (pull out the mains plug).

The compressor in a cabinet requires that the jumper is not mounted in the terminal strip for correct operation (Fig. 24 pos. A). The switch on the cabinet switches the entire product on and off, including the compressor.

If there is a jumper in the terminal strip of the compressor connected to a cabinet, the switch on the compressor cabinet will not work!

If we want to connect a new cabinet to a compressor that previously worked alone (a compressor without a cabinet), then it is necessary to cancel the connection in the compressor terminal strip by performing the procedure, see further. (Fig. 24). (See also chap.13.1. and 13.2).

PROCESS:

Jumper removal (D-C-B-A process)

- Disconnect the product from the mains by pulling the plug out of the socket.
- Remove the cover on the electrical panel (on the compressor).
- The jumper is in the terminal strip – D
- Lift up the fuse cover – C
- Remove the jumper from the terminal strip – B
- Close the fuse cover – A
- Reinstall the cover on the electrical panel.
- Insert the compressor into the box.

- Connect the cabinet to the compressor by inserting the power cord with connector into a power socket. (Fig. 18)
- Connect the device to the electrical mains
- Activate the compressor by turning on the switch on the pressure switch and the switch (4) on the cabinet (Fig. 20).

Compressor with cabinet - when operating in the cabinet, there **must not be** an electrical connection. (Fig. 24 pos. A).

Note: **In case of using a compressor without a cabinet or in case of maintenance, a jumper which is not installed into the terminal strip should be kept** in a designated place next to the electrical panel.

19.10. Cleaning and disinfection of the exterior surfaces of the product

Clean and disinfect the exterior surfaces with neutral cleaning products.

Use of aggressive detergents and disinfectants containing alcohol and chlorides can lead to surface damage and discolouration.

AD dryer maintenance

19.11. Venting pressure from the dryer

The equipment is designed to permit the safe venting of pressure within 10 seconds after the compressor is shut off.

If pressure is not automatically vented from the dryer, then the pressure may be vented manually.

Venting compressed air poses an injury hazard

Wearing hearing protection is recommended given the noise generated by the venting process.

Shut off the compressed air source before venting pressure from the equipment

Venting pressure using the display screen

Press ESC+▼ in the display screen to vent pressure from the equipment.

- First, shut off the compressed air source.
- Then press and hold ESC+▼, which will open all the solenoid valves (inlet and regeneration) for 10 seconds and then vent the pressure from the equipment and connected pneumatic circuits and elements that are not separated from the equipment by the check valve

Manual venting of pressure

- Turn off the compressor.
- Open the vent plugs on the outlet module on the equipment (Fig. 25).


Fig. 25: Venting pressure from the dryer chambers

- Disconnect the hose (2) from the lower part of the condensate separator (3) (Fig. 26).


Fig. 26: Venting pressure from the cooler and condensate separator

The process of manually venting pressure from the equipment is complete after approximately 2 minutes.

19.12. Replacement of the dryer's internal filters

Working with pressurised pneumatic components poses a risk of injury.

Prior to any work, disconnect the equipment from the mains, shut off the compressor and vent all pressure in the equipment to zero.

In normal operation, filter replacement must be performed in the upper part of the dryer at the defined interval.

- Turn off the compressor.
- Check the pressure in the dryer.
- If the dryer chambers are under pressure, proceed in accordance with Chapter 19.11.
- Unscrew the 8 screws (1).
- Disassemble the outlet panel (2) on which the filters (3) are mounted.
- Unscrew the dirty filters (3) and replace with new filters.
- Check the seal (4) on the bottom of the outlet module and replace if necessary.
- Reverse the procedure to reassemble.
- Switch on the compressor.
- Check for any dryer leaks.


Fig. 27: Replacement of internal filters

19.13. Replacement of cassettes with adsorbent media

In normal operation, the replacement of the cassettes with adsorbent media must be performed at the defined interval.

- Turn off the compressor.
- Check the pressure in the dryer.
- If the dryer chambers are under pressure, proceed in accordance with Chapter 19.11.
- Unscrew the 8 screws (1).
- Remove the outlet panel (2).
- Pull out and replace the cassettes (3) with new parts.
- Check the seal (4) on the bottom of the outlet module and replace if necessary.
- Reverse the procedure to reassemble.
- Switch on the compressor.
- Check for any dryer leaks.


Fig. 28: Replacement of cassettes with adsorbent media

19.14. Replacement of the logic valve ball

- Turn off the compressor.
- Check the pressure in the dryer.
- If the dryer chambers are under pressure, proceed in accordance with Chapter 19.11.
- Unscrew the 4 screws (1) and remove the cover (2).
- Remove the ball cover (3).
- Replace the ball (4).
- Check the nozzles (5) and clean as necessary.
- Reverse the procedure to reassemble.
- Check for leaks and the operation of the logic valve and nozzles – check for the cyclical switching of the chambers.


Fig. 29: Replacement of the logic valve ball

19.15. Replacement of the dryer's silencer

Working with pressurised pneumatic components poses a risk of injury.

Operating the equipment without silencers generates high levels of noise. Only replace silencers when the equipment is shut down.

- Unscrew the silencer (1).
- Install a new silencer.


Fig. 30: Replacement of the silencer

19.16. Inspecting the cooled and fan

The equipment, in particular the compressor fan, cooler fan, and the cooler, must be kept clean to ensure efficient drying. Remove dust from the surface of the cooling fins and fans by vacuuming or blowing down with compressed air.

19.17. Replacement of solenoid valves

Risk of electric shock.

Shut off the compressed air source, turn off the equipment and disconnect it from the mains before working on the equipment.

Working with pressurised pneumatic components poses a risk of injury.

Disconnect the equipment from the mains and vent the pressure in the equipment and the pneumatic system to zero before working on the equipment.

In normal operation, the solenoid valves in the bottom of the dryer must be replaced at the defined interval.

- Turn off the compressor.
- Check the pressure in the dryer.
- If the dryer chambers are under pressure, proceed in accordance with Chapter 19.11.
- Unscrew the 1 screw from the valve connector (1).
- Disconnect the valve connector (2).
- Unscrew the 4 screws (3).
- Remove the solenoid valve (4).
- Remove the valve seal (4-1) from the body.
- Physically clean the valve seat surface to remove any impurities.
- Physically clean the 16 screws to remove the thread locking adhesive.
- Install the solenoid valve (Fig. 32).
- Install the new valve seal (4-1).
- Install the new solenoid valve using the 4 screws (3) and use a thread locking adhesive on the threads of the screws (such as Loctite 243).
- Reattach the solenoid valve connector and attach with a screw (1).
- Switch on the compressor.
- Check for any dryer leaks.


Fig. 31: Solenoid valve replacement

Solenoid valve assembly

Replacement solenoid valves are delivered as disassembled replacement parts. The new valve must be assembled before a solenoid valve is replaced.

- Mount the valve coil (4-5) onto the valve body (4-4) and secure with the nut (4-6).
- Insert the valve membrane spring (4-3) into the membrane (4-2) and the insert into the assembled valve coil and body assembly.
- Then mount the valve seal (4-1) onto the dryer body.


Fig. 32: Solenoid valve assembly

19.18. Pressure relief valve

The pressure relief valve automatically begins to vent air from the system if the pressure in the

compressed air circuit exceeds its pre-set value. The pressure relief valve closes as the pressure drops.

The pressure in the compressed air circuit can only increase because of an increase in flow resistance in the compressed air lines or as a result of a dryer malfunction (e.g. solenoid valve malfunction), and therefore the repeated opening of the relief valve requires a dryer function check and repairs if necessary!

Consultation with the manufacturer is required before any adjustment is made to the relief valve!

The outlet openings on the relief valve may not be blocked and the egress of compressed air through them may not be restricted.

2 Dryer pressure valve


Fig. 33: Pressure relief valve

20. LONG-TERM SHUTDOWN

If the compressor is not going to be used for long period, it is recommended to drain all condensate from the air tank and switch on the compressor for about 10 minutes, keeping the drain valve open (1) (Fig. 14). Thereafter switch off the compressor using the switch (2) on the pressure switch (1) (Fig. 20), close the drain valve and disconnect the device from the mains.

AD dryer shutdown

- Vent pressure from the dryer chambers. (Fig. 25).

21. DISPOSAL OF DEVICE

- Disconnect the equipment from the mains.
- Release the air pressure in the pressure tank by opening the drain valve (Fig. 25).
- Dispose of the equipment following all applicable regulations.
- Entrust a specialised company to sort and dispose of waste.
- Worn out components have no negative environmental impact.

TROUBLESHOOTING**Risk of electric shock.****Before interfering with the equipment, first disconnect it from the mains (remove the power socket).****Working with pressurised pneumatic components poses a risk of injury.****Before interfering with the equipment, vent the air tank and the compressed air system to zero pressure.****Troubleshooting may only be performed by a qualified service technician.****Damage to the safety valve could cause pressure to rise to hazardous levels.****Never adjust a safety valve.**

Malfunction	Possible cause	Solution
Compressor does not switch on	No voltage in the pressure switch	Check voltage in socket
		Check circuit breaker switch – switch to position switched-on „I“
		Loose conductor from terminal - repair
		Check electrical cord – replace defective cord
	Motor winding failure, damaged thermal protection	Replace motor or windings
	Capacitor failure	Replace capacitor
Compressor often switches on	Seized piston or other rotating part	Replace damaged components
	Pressure switch does not switch	Check function of pressure switch
	Air leakage in pneumatic system	Check pneumatic system – seal loose connections
	Non-return valve (SV) leakage	Clean non-return valve, replace seals, replace non-return valve
	Large volume of condensed liquid in pressure tank	Drain condensed liquid
Low compressor capacity	Check time of filling the air tank	
Low pressure in the air tank (compressor running constantly)	High air consumption of supplied equipment	Decrease air consumption Use compressor with higher capacity
	Leakage in pneumatic system	Check pneumatic system – seal loose connections
	Low capacity of the pump	Clean / replace the pump
	Pump failure	Clean / replace the pump
	Dryer failure	Replace dryer
Prolonged operation of the	Air leakage in pneumatic system	Check pneumatic system – seal loose connection

compressor	Worn piston ring	Replace worn piston ring
	Dirty inlet filter	Replace dirty filter with new filter
	Incorrect function of solenoid valve	Repair or replace fan or coil
Compressor is noisy (knocking, metal noises)	Damaged piston bearing, piston rods, motor bearing	Replace damaged bearing
	Loose (cracked) dampening element (spring)	Replace damaged spring
Dryer does not dry (condensed water in the air) *	Defective cooler fan	Replace fan Check power source
	Damaged dryer	Replace dryer
	Defective automatic condensate drain	Clean / replace
Degraded drying performance – high-pressure dew point (condensed water in the air)	Low operating pressure	Reduce the demand for air, check the output from the compressed air source, fix any leaks in the distribution system
	Regeneration solenoid valve not working	Check coil operation, replace if damaged Inspect the condition of the valve – clean the valve or replace if problems persist
	Air regeneration nozzle plugged	Clean or replace the nozzle (see product maintenance)
	Cooling fan not working	Check the power source to the fan Replace damaged fan
	Dirty cooler	Inspect the cooler and clean as necessary
	Silencer plugged at outlet from regeneration valve	Inspect the silencers. Clean or replace the silencer if flow resistance is too high or if heavily soiled.
Dryer emitting high levels of noise	Damaged fan	Replace damaged fan
	Damaged silencer	Replace the silencer
	Air leaking through relief valve at dryer inlet	Check the dryer connection to the mains and dryer connections, check the dryer operation, check the dryer's working pressure, and replace defective components.
Air leaking through relief valve at dryer inlet	High working pressure from air source	Check the pressure setting on the compressed air source
	Dryer inlet solenoid valve not working	Check the coil operation, replace if damaged Inspect the condition of the valve – clean the valve or replace if problems persist
	High pressure in equipment resulting from plugged filters	Check the internal filters and accessory filter assemblies. Clean or replace dirty filters.

)* Clean inside surfaces of the air tank thoroughly and remove all condensed liquid after dryer failure.

Once a dryer fault is cleared and after its reassembly, the dryer should quickly be regenerated, best

when using continuous compressor operation at a pressure of around 1 bar for a period of at least 1 hour; then check the dryness of the compressed air.

Check the moisture content of the air exiting the air tank (see the Technical data chapter) to prevent damage to connected downstream equipment.

22. REPAIR SERVICE

Warranty and post-warranty repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

Attention.

The manufacturer reserves the right to make changes to the equipment without notice. Any changes made will not affect the functional properties of the equipment.

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1. Product: (model) DK50 2x4VR/110/M DK50 2x4VR/110S/M		2. Serial number:	
3.1. User's name:			
3.2. Address of installation:			
4. Equipment connected to the compressor:			
5. Installation / Commissioning:		6. Contents of operator training:	
Product completeness check **	Y	Description of the product and functions**	Y
	N		N
Documentation completeness check **	Y	Product operation: turning on/off, controls, control procedures, data on the display panel, alarms, operation in alarm conditions**	Y
	N		N
Installation/connection to equipment **	Y	Product maintenance: maintenance intervals, maintenance procedure, service intervals, operating activities**	Y
	N		N
Functional test **	Y	Safety measures, warnings – their meaning and compliance **	Y
	N		N
Notes:			
7. Operator instructed on safety measures, operations and maintenance:			
Name:		Signature:	
Name:		Signature:	
Name:		Signature:	
8. Installation and instruction performed by: First name/Last name		Signature:	
Company:		Address:	
Phone:			
Email:		Date:	
9. Distributor:			
Company:		Address:	
Contact person:			
Phone:		Email:	

** mark with an "X" in points 5 and 6 (Y - yes /N - no). Enter any observations from points 5 and 6 into the "Notes" section

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